

Diagnostic Automation, Inc.

AccuDiag™

Total T3

ELISA

Test Background

T3 determination is also useful in monitoring both patients under treatment for hyperthyroidism and patients who have discontinued anti-thyroid drug therapy. It is especially valuable in distinguishing between euthyroid and hyperthyroid subjects.

[3144-P1]

The quantitative determination of Total Triiodothyronine concentration in human serum.



Advantages:

- Specific and sensitive technology
- High Throughput Screening (HTS) compatible
- Cost Efficient

Specifications

Principle	Competitive Enzyme Immunoassay
Detection Range	0.15 – 8.0 ng/mL
Sample	50 µL serum/plasma
Incubation Time	75 minutes
Shelf Life	12-14 months from manufacturing date

Sensitivity:

0.04 ng/mL