




Diagnostic Automation/Cortez Diagnostics, Inc.

AccuDiag™





Anti-Phospholipid Screen IgG/IgM ELISA


 www.rapidtest.com
 onestep@rapidtest.com
 818-591-3030

PURPOSE

A test system to screen for the presence of IgG and IgM class autoantibodies against cardiolipin, phosphatidyl serine, phosphatidyl inositol, phosphatidic acid and beta-2-glycoprotein I in human serum or plasma.

ADVANTAGES

-  Very easy to use with little training
-  Highly specific and consistent assay
-  Provides accurate results quickly
-  Reading of results both visually and as absorbance data



SPECIFICATIONS

Anti-Phospholipid Screen IgG/IgM ELISA	
Method	Enzyme Linked Immunosorbent Assay
Principle	Indirect; Antigen Coated Plate
Detection Range	Semi-Quantitative Positive, Negative Control and Cut-off
Sample	10 µL Serum
Total Time	~ 65 min.
Shelf Life	12 Months from the manufacturing date

SENSITIVITY

IgG: 91.8%
IgM: 53.4%

SPECIFICITY

IgG: 97.3%
IgM: 96.7%

Contact us today for more information!