




Diagnostic Automation/Cortez Diagnostics, Inc.

AccuDiag™

TBG

(Thyroid Binding Globulin)

ELISA

 www.rapidtest.com
 onestep@rapidtest.com
 818-591-3030

PURPOSE

Intended for quantitative determination of Thyroxine Binding Globulin (TBG) concentration in human serum, plasma or whole blood.

ADVANTAGES

-  Very easy to use with little training
-  Highly specific and consistent assay
-  Provides accurate results quickly
-  Reading of results both visually and as absorbance data

SPECIFICATIONS

TBG (Thyroid Binding Globulin) ELISA	
Method	Enzyme Linked Immunosorbent Assay
Principle	Competitive Immunoassay
Detection Range	Quantitative
Sample	10 µL Serum/Plasma/WB
Total Time	~ 45 min.
Shelf Life	12 Months from the manufacturing date

SENSITIVITY
1.0 µg/mL

Contact us today for more information!