

AccuDiag™ CA-19-9 Cancer ELISA

[6909-P1]

Intended Use

The CA-19-9 ELISA Test is intended for the quantitative determination of CA-19-9 in human serum.

A group of mucin-type glycoprotein Sialosyl Lewis Antigen (SLA), such as CA-19-9 and CA-19-5, have come to be recognized as circulating cancer associated antigens for gastrointestinal cancer.



Advantages

- Sensitive and specific technology
- Runs up to 96 tests
- Great for High Throughput Screening (HTS)
- Results both visually and as absorbance data
- Easy procedural steps

Specifications

Principle	Immunoenzymometric Sequential Assay
Detection Range	0–500 U/mL
Sample	25 µl serum
Sensitivity	1.0 U/mL
Incubation Time	135 minutes
Shelf Life	12 Months from the manufacturing date